

Linak Battery Charger for Lifter



- Functional and elegant design
- Approx. charging time for BAJ1/BAJL battery pack: 5 hours
- The charger is of the SMPS type (Switch Mode Power Supply)
- Mains voltage: 100-240 V / 50-60 Hz
- Yellow light diode indicates whether the battery is charging
- Greenlight diode indicates whether the charger is connected to the mains
- Charging current up to 650 mA
- Protection class: IPX5 (only when mounted vertically)
- Colour: Bottom dark grey (RAL 7016) and top light grey (RAL 7035)
- 3.2 m exchangeable mains cable (black straight)

Code	Display Name
NZ-MBT401210	Linak Battery Charger for Lifter