

## Aspire Commode Open Seat



- Moulded and seamless soft PU seats designed for comfort and durability
- Easy clean, fully sealed surface for infection control
- Choose the size appropriately to fit your Aspire Shower Commode
- Adaptable to most shower commode brands
- M6 Tri Nut assembly kit needed for seat attachment. Comes as standard with all commodes. Available as a spare part if lost (GSP342405)

### Warranty Information

- 1 Year



Code	Display Name
NZ-BTC066600	Open Front Commode Seat - 460 mm Wide
NZ-BTC066750	Open Front Commode Seat - 530 mm Wide
NZ-BTC066800	Open Front Commode Seat - 600 mm Wide