

Aspire STAR Mover



- The Aspire Star Mover offers a safe and comfortable solution for patients and carers doing sit to stand transfers
- The compact and versatile design is perfect for Homecare, Institutional and Aged Care environments
- Wide seat entry suits a larger range of body shapes
- Low footplate height makes patient entry and exit easier
- Large 100 mm rear locking castors have low rolling resistance and offer a smooth ride during transportation
- Narrow base improves access in tight spaces
- Height adjustable knee support ensures correct positioning for a larger patient height range
- Easy access seat pad handles provide reach-around access to carers
- Soft touch grab bar provides a comfortable patient support for improved stability and confidence during transfers
- Width adjustable legs with tapered close position allow access between and/or around chair legs
- Knee pad height: 310 - 580 mm
- Seat pad height: 850 mm
- Handle height: 1025 mm
- Leg width range: 520 - 740 mm
- Footplate height: 70 mm
- Turning Circle: 1055 mm

Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
NZ-LSS390580	Standing & Transport Aid - Aspire STAR Mover - Adjustable Legs - 185kg	920 mm	615 - 920 mm	1030 mm	185 kg