

Aspire Lifecomfort Deluxe Mattress



The LIFECOMFORT® DELUXE includes three generous foam layers; visco-elastic memory foam upper; high resilience comfort layer; and a supportive base. An ideal solution for Aged Care it provides users with an immersive, pressure relieving support surface for lasting comfort. Fitted with superior, hospital-grade polyurethane coated covers designed to withstand the rigors of usage and cleaning requirements in Home Care and Aged Care

- **Durable Hospital Grade Foam** that is treated for antimicrobial and antibacterial properties
- **Breathable, 4-Way Stretch Waterproof Cover** to reduce shear, friction and maintain skin microclimate
- **Medical Grade Waterfall Cover** is fully-welded over the zipper to maximise infection control
- **Asset Management Ready** with unique batch and serial numbers printed on all covers
- **Carry Handle** - Integrated handles for safe manual handling
- **Hospital Grade** - Suitable for bleach concentrations up to 10,000 ppm
- **4-sided Zipper** - Easily removed for effective cleaning
- **4-way Stretch** - For shear mitigation and friction reduction
- Immersive surface contours and cradles around the user for personalised pressure relief and comfort
- Machine Washable Covers: Up to 95 degrees
- Engineered to move with adjustable beds
- **Foam Layers:**
 - IMMERSIVE SURFACE - Immersive visco-elastic
 - COMFORT LAYER - High resilience foam for support and pressure redistribution
 - SAFETY & SUPPORT BASE - Strong base support for maximum durability

- **Applications** - Pressure Care: Medium
 - Aged Care
 - Ambulant patients
 - Home Care

Warranty Information

- 5 Years Foam
- 4 Years Cover

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
NZ-BEM046510	Lifecomfort Deluxe Mattress - Single	2000 mm	900 mm	150 mm	30 - 150 kg	Low Risk
NZ-BEM046455	Lifecomfort Deluxe Mattress - King Single	2030 mm	1050 mm	150 mm	30 - 150 kg	Low Risk