

Aspire Deluxe Hygiene Access Sling with Head Support - Clip - Polyester



- Clip fitting option to fit sling to a pivot frame
- Australian standards tested and registered with the TGA
- This sling offers an adjustable head support for additional patient safety
- Intuitively designed, this sling offers generous underarm padding, quilted leg straps and a reinforced wide-cut back for comfort. The chest strap features velcro and buckle fastening for easy positioning and added safety.
- Padded poly fabric & generous under arm padding for increased patient comfort
- Versatile sling offering body contoured fabric shaping for maximum patient comfort and support
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable - loop webbing has a breaking strain in excess of one tonne and is designed so that it will not tangle or “cork-screw” during cleaning
- User Safety - smarts feature such as velcro/buckle combinations for simplified positioning and greater user safety
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- For FIM score 2 or 3

Warranty Information

- 1 Year

Code	Display Name	Safe Working Load	Material
NZ-LSS390955	Deluxe Hygiene Access Sling with Head Support - Clips - Polyester - Extra Large	300 kg	Polyester