

## EasyMove Bed Transport System



- The EasyMove Transport System is an ideal transport solution designed to minimise manual handling.
- It allows effortless lifting and tilting of the bed to its transport position, while the large 100mm castors offer smooth and easy manoeuvring for one person.
- The EasyMove is also a perfect maintenance platform facilitating easy access to mechanical parts under the bed platform.

### Warranty Information

- 1 Year



Code	Display Name
NZ-BEA030055	Bed Mover/Transporter - Aidacare FLX & FL250 Bed (Set of 2))