

## Aspire Deluxe Standing Sling - Polyester - Black Clips



- Clip fitting option to fit sling to a pivot frame
- Australian standards tested and registered with the TGA
- Intuitively designed, this sling offers generous underarm padding and a reinforced wide-cut back for comfort
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable keyhole clip attachment points for safety and security during transfers
- User Safety - smarts feature such as velcro/buckle combinations for simplified positioning and greater user safety
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- For FIM score 3 or 4

### Warranty Information

- 1 Year

| Code          | Display Name   | Safe Working Load | Material  |
|---------------|--|-------------------|-----------|
| NZ-LSS390920B | Deluxe Padded Standing Sling - Polyester - Black Clips - Small | 200 kg            | Polyester |

| Code          | Display Name   | Safe Working Load | Material  |
|---------------|--|-------------------|-----------|
| NZ-LSS390925B | Deluxe Padded Standing Sling - Polyester - Black Clips - Medium      | 200 kg            | Polyester |
| NZ-LSS390930B | Deluxe Padded Standing Sling - Polyester - Black Clips - Large       | 300 kg            | Polyester |
| NZ-LSS390935B | Deluxe Padded Standing Sling - Polyester - Black Clips - Extra Large | 300 kg            | Polyester |