


Aspire Deluxe Standing Sling - Polyester - Black Clips



- Clip fitting option to fit sling to a pivot frame
- Australian standards tested and registered with the TGA
- Intuitively designed, this sling offers generous underarm padding and a reinforced wide-cut back for comfort
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable keyhole clip attachment points for safety and security during transfers
- User Safety - smarts feature such as velcro/buckle combinations for simplified positioning and greater user safety
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- For FIM score 3 or 4

Warranty Information

- 1 Year

Code	Display Name	Safe Working Load	Material
NZ-LSS390920B	Deluxe Padded Standing Sling - Polyester - Black Clips - Small	200 kg	Polyester

Code	Display Name	Safe Working Load	Material
NZ-LSS390925B	Deluxe Padded Standing Sling - Polyester - Black Clips - Medium	200 kg	Polyester
NZ-LSS390930B	Deluxe Padded Standing Sling - Polyester - Black Clips - Large	300 kg	Polyester
NZ-LSS390935B	Deluxe Padded Standing Sling - Polyester - Black Clips - Extra Large	300 kg	Polyester