

Aidacare Aspire 200S Steel Standing Lifter



The Aidacare Aspire 200S is a robust and durable steel standing lifter capable in the most demanding facility environments

- **Bullhorn Lifting** Handles for patient stability and 6 sling attachment points
- **Ergonomic Handles** Provide multiple grip points to suit carers of different heights and strengths
- **Ergonomic Handset** The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- **Powered Leg Function** Opens and closes the legs for lifting stability or accessibility
- **Sturdy steel frame** with industry leading Linak electronics and warranty provision
- **Safety** Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- **Adjustable Padded Legrests** Provide comfort and support during transfers & can be adjusted to suit the patient
- **Intelligent Diagnostics** State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- **Tente Casters:** 2 x 75mm front castors for maximum under bed clearance and accessibility. 2 x 75mm rear locking castors for safety and manoeuvrability
- Item Weight: 44.6 kg
- Tested and compliant to AS/NZS ISO 10535:2011 (the Australian standard for hoists)

Warranty Information

- Frame: 5 Years
- Electronics: 3 Years

- Battery and Handset: 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range	Electric Leg
NZ-LSS390650	Aidacare Aspire 200S Standing Patient Lifter	1150 mm	656 - 1016 mm	1790 mm	200 kg	Steel	44.6 kg	850 - 1740 mm	575 - 936 mm	True