

Aspire Overbed Table - Laminate Flat Top



- Lightweight mobile table positions overbed or chair
- Can be used for eating, drawing or other activities
- Flat top ideal for hospital, facility or home use
- Full height adjustable with a fixing knob for additional security
- C base shape fitted with four castors enables close positioning to the bedside or chair
- Powder coated mild steel frame construction
- Cream standard laminate top with vinyl wrap around the edge, or thermoform cream top which is heat resistant
- Available in White & Beech Laminate coloured tops

Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Total Height
NZ-BEA006701	Aspire Overbed Table - Laminate Flat Top - Beech	430 mm	800 mm	730 - 1095 mm