

## Active Air 8 Alternating Mattress Replacement



The Active Air 8 is a dynamic mattress replacement suitable for users at “High Risk” of Pressure Injury. Featuring AUTODETECT, and AUTOPROTECT, together with CELL-IN-CELL technology, the Active Air 8 provides the highest levels of care into Acute and Long Term care environments. The Active Air 8’s innovative features ensure lowest interface pressures with an easy-to-use digital pump.

- Digital Pump – The Easy-to-read LCD screen provides clear and precise information on cycle time, operating mode and fault codes to assist with troubleshooting
- AutoLock – Control panel automatically locks after 30 seconds, preventing accidental mode changes
- AutoDetect – Based on the user’s weight and size, appropriate pressure output is automatically set, ensuring optimal care and protection
- AutoProtect – When the head of the bed is raised to 45° or more, air pressure automatically increases to prevent bottoming out
- Continuous Low Pressure (CLP) – When activated, this mode increases user immersion and envelopment
- 24hr Transport Mode – In the event of a power failure, this feature provides a security buffer allowing patients to remain on an inflated surface until power is restored
- Sleep Quality - Innovative pump technology reduces noise and vibration, for a quieter sleep environment
- Infection Control - Welded seams and zip protection prevent fluid and bacterial ingress into the system
- Skin Microclimate - The multi-stretch cover reduces shear and offers a high MoistureVapour transfer rate (MVTR)
- Workplace Health & Safety - The power cord is retained within the overlay eliminating its potential as a trip hazard

- Pressure Care Outcomes - Intelligent cell design maximises the duration of pressure reduction in the alternation cycle
- Laundry & Hygiene - Machine washable cover at 95° to achieve superior sanitising

## Warranty Information

- Pump: 2 Years
- Replacement: 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
NZ-PTM639300	Active Air 8 Alternating Mattress Replacement - Single	2000 mm	900 mm	200 mm	30 - 200 kg	High Risk
NZ-PTM639310	Active Air 8 Alternating Mattress Replacement - King Single	2000 mm	1050 mm	200 mm	30 - 200 kg	High Risk