

Aspire Assist 3 Manual Wheelchair



- **The Assist 3 Wheelchair** – is designed for an individual or facility who requires a robust chair that will meet the rigorous demands of the community and healthcare environments
- **Swing away and removable arm rests** – this facilitates user transfers and also makes for easier transport and storage
- **Height adjustable and removable foot rests** – this facilitates user transfers, as well as allowing for optimal positioning for the changing needs of the user
- **Heel loops** on the foot rests – increases user safety
- **Quick release 24" solid rear tyres** – the wheels are easily removed to reduce the weight of the chair for transport and storage
- **An easy to use, rear wheel brake** with a large handle for ease of access

This wheelchair is NOT crash-tested and is NOT certified for transport use in vehicles.



Warranty Information

- 2 Years Frame
- 1 Year Components

Code	Display Name	Safe Working Load	Material	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight	Backrest Height
NZ-MWS452200	Aspire Assist3 Folding Wheelchair - 350mm	140 kg	Steel	Silver	400 mm	350 mm	460 - 500 mm	460 - 500 mm	17.3 kg	440 mm
NZ-MWS452210	Aspire Assist3 Folding Wheelchair - 400mm	140 kg	Steel	Silver	400 mm	400 mm	460 - 500 mm	460 - 500	17.5 kg	440 mm

Code	Display Name	Safe Working Load	Material	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight	Backrest Height
NZ-MWS452220	WChair Folding - Aspire Assist3 (Self Propelled) - 450mm - 140kg - Silver	140 kg	Steel	Silver	400 mm	450 mm	460 - 500 mm	460 - 500 mm	17.7 kg	440 mm
NZ-MWS452230	WChair Folding - Aspire Assist3 (Self Propelled) - 500mm - 140kg - Silver	140 kg	Steel	Silver	400 mm	500 mm	460 - 500 mm	460 - 500 mm	17.9 kg	440 mm