

## Aspire Low Back Classic Day Chair



- A versatile, low back orthopaedic chair
- Contoured Backrest – Provides postural support
- Softly padded armrests, seat and back upholstery are comfortable and supportive
- Non-Slip Rubber Tips – Provide stability whilst sitting or transferring
- Easy clean vinyl surfaces
- Height Adjustable – Easy push button adjustment for comfortable positioning and changing needs
- Lightweight Design - lightweight yet sturdy design for easy manoeuvrability
- Powder Coated Steel Frame – Easy to clean, strong and durable
- Available in either Slate or Champagne Vinyl

### Warranty Information



- 1 Year

Code	Display Name	Safe Working Load	Material	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
NZ-CHP208000	Aspire Low Back Classic Day Chair - Champagne Vinyl	160 kg	Steel, foam, vinyl	Champagne Vinyl	440 mm	460 mm	430 - 605 mm	430 - 605 mm	440 mm
NZ-CHP208005	Aspire Low Back Classic Day Chair - Slate Vinyl	160 kg	Steel, foam, vinyl	Slate Vinyl	440 mm	460 mm	430 - 605 mm	430 - 605 mm	440 mm