

Aspire Over Toilet Aid - MAXI Adjustable



- Super heavy duty adjustable width over toilet aid
- Soft padded flip up seat for easy cleaning
- Reinforced aluminium frame accommodates bariatric users and is corrosion resistant
- Height adjustable legs with easy circlip adjustment points
- Raised armrests to assist the user in transferring on and off the stool
- Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability
- To create a MAXI bedside commode (BEC045120), simply add a padded backrest (as shown in MAXI Shower Chair above)



Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
NZ-BTT145020	Aspire MAXI Over Toilet Aid	645 - 845 mm	670 - 790 mm	310 kg	455 mm	575 - 775 mm	455 - 575mm	455 - 575mm