

Revolution Chair



- The Revolution Chair incorporates a unique, patented “swivel-slide” mechanism.
- Swivel – to allow easy access to the chair seat.
- Slide – to bring you close to the table edge.
- At all times, the chair legs remain fixed and stable.
- Sturdy steel frame, waterfall polyurethane arms.
- Bleach-cleanable commercial grade vinyl seat and back.

Code	Display Name	Total Width	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
NZ-CHP196832	Revolution Chair - Height Adjustable - Onyx	505 mm	120 kg	Onyx	485 mm	450 mm	450 - 530 mm	450 - 530 mm
NZ-CHP196834	Revolution Chair - Height Adjustable - Sand	505 mm	120 kg	Sand	485 mm	450 mm	450 - 530 mm	450 - 530 mm
NZ-CHP196830	Revolution Chair - Height Adjustable (Custom - Specify Colour)	505 mm	120 kg	-	485 mm	450 mm	450 - 530 mm	450 - 530 mm
NZ-CHP196840	Revolution Chair - Bariatric (Custom - Specify Colour)	800 mm	250 kg	-	450 mm	660 mm	480 mm	480 mm