

## Aidacare Aspire A200S Aluminium Standing Lifter



The Aidacare Aspire A200S is a compact, lightweight and manoeuvrable standing lifter ideal for home, community and facility environments.

- **Adjustable Depth Yoke:** Yoke is depth adjustable for optimal patient fit
- **Ergonomic Handles:** Provide multiple grip points to suit carers of different heights and strengths
- **Ergonomic Handset:** The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- **Powered Leg Function:** Opens and closes the legs for lifting stability or accessibility
- **Build Quality:** Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision
- **Safety:** Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- **Adjustable Padded Legrest:** Provide comfort and support during transfers & can be adjusted to suit the user
- **Intelligent Diagnostics:** State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- **Tente Casters:** 2 x 75mm front castors for maximum under bed clearance and accessibility. 2 x 100mm rear locking castors for safety and manoeuvrability.
- Item Weight: 30.4 kg
- Tested and compliant to AS/NZS ISO 10535:2011 (the Australian standard for hoists)

### Warranty Information

- Frame: 5 Years
- Electronics: 3 Years

- Battery and Handset: 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range	Electric Leg
NZ-LSS390640	Aidacare Aspire A200S Standing Patient Lifter	985 mm	590 - 915 mm	1455 mm	200 kg	Aluminium	30.4 kg	890 - 1395 mm	490 - 815 mm	True