

Aspire Bed Cradles



- The Aspire bed cradle is designed to keep bed linen off sensitive areas of the patient's lower leg or feet by creating clearance between bed linen and the mattress
- The lower frame slides underneath the mattress, securing the cradle to the bed, while the upper frame is positioned horizontally over the lower limbs
- Bedding is then placed over the frame to remove any weight or pressure from the sensitive limbs below
- Height adjustable frame allows the user to customise the space between mattress surface and bedding (BEA014610 only)
- Powder Coated Mild Steel Construction

Warranty Information

- 1 Year

Code	Display Name	Total Height	Material
NZ-BEA014610	Aspire Bed Cradle - Height Adjustable	440 - 640 mm	Steel