

Aspire Tilt & Recline Shower Commode



- 30° of backrest angle adjustment in 10° increments for effortless repositioning and maximum user comfort
- 0° to 35° of gas-assisted tilt-in-space functionality minimises manual handling efforts during the showering process
- 4x Tente total locking castors provide brake access from all angles
- Open back straps provide greater access for bathing needs
- Backrest may be configured as full length or 3/4 length to accommodate an advanced headrest
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Height adjustable, swing away footrests with calf pads for optimal positioning, support and ease of transfers
- Tool-free height adjustment to 4 settings allow the commode to be easily adapted to varying needs
- Overall Height: 1315 - 1485 mm
- Overall Width: 560 mm (BTC066030), 530 mm (BTC066040)
- Seat to Backrest: 900 mm
- Anterior Seat Height: 460 - 630 mm
- Posterior Seat Height: 460 - 620 mm
- ***Seat sold separately**
- Note: Seat is supplied automatically with any rental equipment

Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

Code	Display Name	Safe Working Load	Seat Width	Backrest Height
NZ-BTC066030	Tilt & Recline Commode - Attendant Propelled - 460mm Wide	200 kg	460 mm	900 mm
NZ-BTC066040	Tilt & Recline Commode - Attendant Propelled - 530mm Wide	200 kg	530 mm	900 mm