

## Aspire Lifecomfort Evacuation Sheet



- Designed for the emergency evacuation of the non-ambulant person, the Aidacare Lifecomfort Evacuation Sheet
- Fits discreetly underneath the resident's mattress
- Allows bed making to occur as normal
- Remains readily accessible in case of emergency evacuation
- Is fire retardant and highly durable
- The Evacuation Sheet was designed to help rescue teams, carers or nurses to slide the mattress and resident along any floor surface or even downstairs, increasing resident safety and preventing injury to nursing staff

### Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Safe Working Load
NZ-TSE672526	Lifecomfort Evacuation Sheet	1900 mm	900 mm	300 kg