

Aspire COVE Pressure Relief Chair



- The Aspire Cover chair helps to improve the postural, pressure care and comfort outcomes for users that require medium to high level support.
- Users can be effortlessly moved between rooms at home, or in care environments. This offers improved inclusivity and social interaction.
- **Anatomical Legrest Movement:** As the legrest raises, it extends outwards tracking knee joint range of motion and maximises lower limb support. When combined with correct seat depth adjustment, taller users maintain access under tables promoting social inclusiveness and dining independence.
- **Infinite Tilt Adjustment:** The gas assisted tilting mechanism provides a smooth and intuitive motion to maximise user comfort and confidence. By tilting up to 40 degrees, sitting pressure is redistributed from vulnerable bony prominences to the lower back and assists in maintaining postural alignment.
- **Configurable Pressure Support Cushioning:** The integrated seat cushion combines layered memory foam with a segmented hybrid air cell that offers convenient pressure adjustment to carers and maintenance staff. Reinflation is as simple as opening and closing a valve.
- **Ease of Movement:** High performance TENTE Levina castors have been chosen to reduce push/pull forces during transportation. This improves workplace health and safety outcomes for care staff and family members.
- Directional locking castors at the front of the chair assist during transit, particularly on uneven flooring or carpets with pronounced grains.
- **Infection Control:** The sophisticated medical grade upholstery has been chosen to perform under modern infection control protocols. Removable seating elements and smooth wipeable surfaces simplify cleaning processes between users. The biostatic patient surface prevents bacteria and virus ingress.
- **Market Leading Upholstery:** The patient surface is genuine DARTEX made in the UK. This breathable highly elastic surface maximises skin microclimate and pressure redistribution.

- External surfaces are medical grade coated fabrics from Warwick Australia.
- Hardy yet stylish patterns offer an appealing aesthetic and blend seamlessly in with modern Aged Care interior designs.
- The layered seat structure accommodates after-market pressure care products such as ROHO cushions, without compromising chair performance or infection control outcomes (purchased separately).
- Transitional pressure management armrests.
- Ergonomic carer push handles.
- Tool free adjustments.
- Aspire Cove Chair Accessories available here. {insert link to CHP137}

Warranty Information

- Fabric: 1 year
- Frame: 3 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
NZ-CHP198978	Chair Pressure Relief - Aspire COVE Chair - Teal Vinyl - Medium - 160kg	1640 mm	780 mm	1270 mm	160 kg	Teal Vinyl	500 - 580 mm	550 mm	540 mm	540 mm	740 mm
NZ-CHP198979	Chair Pressure Relief - Aspire COVE Chair - Teal Vinyl - Large - 160kg	1690 mm	830 mm	1300 mm	160 kg	Teal Vinyl	500 - 580 mm	600 mm	540 mm	540 mm	740 mm