

Alerta - Wireless/AlertaDetect - Alarm Monitor



- The Alerta Wall Point Receiver can be connected to a nurse call system for centralised activation, or used on its own, independent of a nurse call system to give localised alerts through this device only.
- Volume can be deactivated for silent in-room operation when connected to a nurse call system.
- Can be wirelessly paired with up to 5 Alerta wireless products.
- Compatible with the Deluxe Alertamat, Floor+ Alertamat, Bed Alertamat, Chair Alertamat, Alerta Detect and the Alerta Button.
- It uses radio frequency technology to connect, with a wireless range of up to 95 metres.
- Once paired with Alerta wireless devices it only transmits to those devices and will not interfere with other Alerta Wall Point Receivers or devices in use nearby.
- Wireless connect-and-play system.
- No cables, extra versatile.
- High, low and silent volume control.
- Low battery alert.
- Mains or battery powered.
- Takes 2x 1.5v AA batteries.

Code	Display Name
NZ-BEA007670	Alerta - Wireless/AlertaDetect - Alarm Monitor