

## Aspire Lifecomfort Pressure Reducing Mattress



The LIFECOMFORT® PRESSURE REDUCING does exactly that. With visco-elastic memory foam upper, plush comfort and resilient support layers and a strong base it provides four tiers of immersive comfort and support for users requiring additional pressure relief. Integrated sidewalls create roll resistance at the bed edge and provide extra rigidity to ease bed entry and exit. Fitted with superior, hospital-grade polyurethane coated covers designed to withstand the rigors of usage and cleaning requirements in Home Care and Aged Care.

- **Durable Hospital Grade Foam** that is treated for antimicrobial and antibacterial properties
- **Breathable, 4-Way Stretch Waterproof Cover** to reduce shear, friction and maintain skin microclimate
- **Medical Grade Waterfall Cover** is fully-welded over the zipper to maximise infection control
- **Asset Management Ready** with unique batch and serial numbers printed on all covers
- **Carry Handle** Integrated handles for safe manual handling
- **Hospital Grade** - Suitable for bleach concentrations up to 10,000 ppm
- **4-sided Zipper** - Easily removed for effective cleaning
- **4-way Stretch** - For shear mitigation and friction reduction
- Immersive surface contours and cradles around the user for personalised pressure relief and comfort
- Machine Washable Covers: Up to 95 degrees
- \*Maximum user weight 300 kg (150kg per user).
- **Foam Layers:**
  - IMMERSIVE SURFACE - Immersive visco-elastic
  - COMFORT LAYER - Plush support layer catering to smaller, lighter weight users

- SUPPORT LAYER - High resilience foam for support and pressure redistribution
- SAFETY & SUPPORT BASE - Strong base support for maximum durability
- INTEGRATED SIDE WALLS - Strong rigid side layers for transfer efficiency and safety
- **Applications** - Pressure Care: Medium
  - Aged Care
  - Ambulant patients
  - Home Care
  - Palliative Care

## Warranty Information

- 5 Years Foam
- 4 Years Cover

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
NZ-BEM046610	Lifecomfort Pressure Reducing Mattress - Single	2000 mm	900 mm	150 mm	30 - 150 kg	Medium Risk
NZ-BEM046660	Lifecomfort Pressure Reducing Mattress - King Single	2030 mm	1050 mm	150 mm	30 - 150 kg	Medium Risk