

Aidacare Medication Cart



- The Aidacare Medication Cart is the powerful all-round solution for a multitude of applications including Treatment, Procedure, IV, Trauma, Dressing, Supply and General cart
- Aidacare carts are designed to offer an attractive, practical and durable cart with multiple configurations to enable the efficient organisation and delivery of medication and equipment
- Robust panels to minimise dents and scratches
- Large work area with slide out work surface are standard
- Easy roll 125 mm castors
- Fully customisable with a host of accessories and drawer configurations.
- **Standard Inclusions:** Drawer Dividers and Extended Sliding Surface
- Lock System Options:
 - **Key lock:** The key lock system is a proven method of securing your medical cart and providing the level of security necessary
 - **Digital lock:** A user friendly system where the cart is locked and unlocked using a pre-set code ensuring ease of access and the security required in clinical environments. The digital lock is robust and waterproof
- Standard drawer configuration detailed below
- **Small Medication Cart:**
 - 1 small 3' drawer (7 cm)
 - 2 medium 6' drawers (16 cm)
 - 1 large 9' drawer (25 cm)
- **Medium Medication Cart:**
 - 2 small 3' drawers (7 cm)
 - 2 medium 6' drawers (16 cm)
 - 1 large 9' drawer (25 cm)
- **Large Medication Cart:**
 - 1 small 3' drawer (7 cm)

- 3 medium 6' drawers (16 cm)
- 1 large 9' drawer (25 cm)

Warranty Information

- Structural: 5 Years Warranty
- Digital Lock: 2 Years Warranty
- Accessories: 1 Year Warranty

Code	Display Name	Total Length	Total Width	Total Height	Colour
NZ-WAT909510	Aidacare Medication Cart - Small - Digital Lock - Cream	790 mm	610 mm	900 mm	Cream
NZ-WAT909560	Aidacare Medication Cart - Small - Key Lock - Cream	790 mm	610 mm	900 mm	Cream
NZ-WAT909540	Aidacare Medication Cart - Medium - Digital Lock - Cream	790 mm	610 mm	1000 mm	Cream
NZ-WAT909590	Aidacare Medication Cart - Medium - Key Lock - Cream	790 mm	610 mm	1000 mm	Cream
NZ-WAT909460	Aidacare Medication Cart - Large - Digital Lock - Blue	790 mm	610 mm	1100 mm	Blue
NZ-WAT909550	Aidacare Medication Cart - Large - Digital Lock - Cream	790 mm	610 mm	1100 mm	Cream