

Aspire Overbed Table - Gunmetal Frame



- Lightweight mobile table positions overbed or chair
- Can be used for eating, drawing or other activities
- Flat top ideal for hospital, facility or home use
- Full height adjustable with a fixing knob for additional security
- C base shape fitted with four castors enables close positioning to the bedside or chair
- Powder coated mild steel frame construction in Gunmetal Grey
- Tabletop 800(W) x 430(L) mm
- Internal base 670 mm



Code	Display Name	Total Height	Colour
NZ-BEA006705	Overbed Table - Gunmetal Frame - Laminate Top - Sublime Teak	730 - 1095 mm	Sublime Teak
NZ-BEA006706	Overbed Table - Gunmetal Frame - Laminate Top - Milano Walnut	730 - 1095 mm	Milano Walnut