

Alerta - Wireless - Detect Sensor System



no image available

- Infrared beam senses when a patient leaves their bed, alerting carers to minimise potential falls risk
- Can be used as a standalone system (sensor & receiver) to create an alarm within the room or directly placed into nurse call system (sensor only)
- Wireless plug and play ensures carer can quickly set up and operate with ease
- System includes: Sensor & Receiver, Stereo & Mono Nurse Call Cables
- Wireless solution that removes excess cables and provides a streamline falls prevention solution
- Large sensor range ensures all types of room configurations are covered
- Mains power cable is available as an optional accessory & batteries are not included

Warranty Information

- 1 Year



Code	Display Name
NZ-BEA007610	Alerta - Wireless - Detect Sensor System (Includes Sensor, Receiver)