


Aspire Deluxe Standing Sling - Polyester



- Australian standards tested and registered with the TGA
- Intuitively designed, this sling offers generous underarm padding and a reinforced wide-cut back for comfort
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable - loop webbing has a breaking strain in excess of one tonne and is designed so that it will not tangle or “cork-screw” during cleaning
- User Safety - smarts feature such as velcro/buckle combinations for simplified positioning and greater user safety
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°c temperature
- For FIM score 3 or 4

Warranty Information

- 1 Year

Code	Display Name	Safe Working Load	Material
NZ-LSS390720	Lifter Sling - Aspire - SML - Deluxe Padded Standing - Poly - 200kg	200 kg	Polyester

Code	Display Name	Safe Working Load	Material
NZ-LSS390725	Lifter Sling - Aspire - MED - Deluxe Padded Standing - Poly - 200kg	200 kg	Polyester
NZ-LSS390730	Lifter Sling - Aspire - LGE - Deluxe Padded Standing - Poly - 300kg	300 kg	Polyester
NZ-LSS390735	Lifter Sling - Aspire - XLGE - Deluxe Padded Standing - Poly - 300kg	300 kg	Polyester
NZ-LSS390736	Deluxe Standing Sling - Polyester - XXL	300 kg	Polyester (Loops)
NZ-LSS390737	Deluxe Standing Sling - Polyester - XXXL	300 kg	Polyester (Loops)