

## Active Air 5 Alternating Overlay



The Active Air 5 provides an effective solution for higher weight users considered to be “At Risk” of Pressure Injury. Featuring AUTO DETECT and INTELLIGENT STATIC MODE which provide lower interface pressure and enhanced comfort, the Active Air 5 is suitable for use in a wide range of care settings from Home Care and Acute to Long Term care.

- Digital Pump – The Easy-to-read LCD screen provides clear and precise information on cycle time, operating mode and fault codes to assist with troubleshooting
- AutoLock – Control panel automatically locks after 30 seconds, preventing accidental mode changes
- AutoDetect – Based on the user’s weight and size, appropriate pressure output is automatically set, ensuring optimal care and protection
- Continuous Low Pressure (CLP) – When activated, this mode increases user immersion and envelopment
- 24hr Transport Mode – In the event of a power failure, this feature provides a security buffer allowing patients to remain on an inflated surface until power is restored
- Sleep Quality - Innovative pump technology reduces noise and vibration, for a quieter sleep environment
- Infection Control - Welded seams and zip protection prevent fluid and bacterial ingress into the system
- Skin Microclimate - The multi-stretch cover reduces shear and offers a high MoistureVapour transfer rate (MVTR)
- Workplace Health & Safety - The power cord is retained within the overlay eliminating it's potential as a trip hazard
- Pressure Care Outcomes - Intelligent cell design maximises the duration of pressure reduction in the alternation cycle
- Laundry & Hygiene - Machine washable cover at 95° to achieve superior



sanitising

## Warranty Information

- Pump: 2 Years
- Overlay: 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
NZ-PTM639200	Active Air 5 Alternating Overlay - Single	2000 mm	850 mm	130 mm	30 - 175 kg	Medium Risk
NZ-PTM639210	Active Air 5 Alternating Overlay - King Single	2000 mm	1050 mm	130 mm	30 - 175 kg	Medium Risk
NZ-GSP493030	Replacement Pump - Aspire Active Air 5 Mattress Overlay - Single	-	-	-	-	-