

## Aidacare FLX Floorline Bed



- Ultra Low Height: Market leading 100mm low height whilst maintaining full ground clearance
- Attendant Control Panel: Allows caregivers functional control and the ability to lock out certain bed functions from being operated via the handset
- Underbed Clearance: Full clearance under the bed minimises the risk of cord impingement and protects flooring from damage
- Custom Height Function: Set the perfect height for the individual user. Caters to a large range of user heights
- SmartDescent: Bed lowering speed is reduced once the bed reaches 300mm maximizing carer potential to remove obstacles as the bed lowers to the floor
- Underbed Lighting: Creates a gentle halo around the bed to reduce falls risk and improve carer safety
- Backrest Auto-Regression: Reduces shear forces on the spine and sacral regions. Minimizes abdominal compression when back is raised supporting cardio pulmonary and respiratory function
- One Touch Chair Positioning: Improves social interaction, assists respiratory and cardio pulmonary function and creates a comfortable position for sitting, dining and transfers
- Integrated Bed Extension: 100-200mm length extension is easily operated and ensures a large range of user heights will fit comfortably without compromise
- Trendelenburg (and Reverse): Augments patient positioning, comfort and assists with bed mobility\*.
- Quick Swap Actuators: Allow for quick and easy inspection and maintenance, minimising downtime
- TiMotion Electronics: Medical grade electronics from a leading global manufacturer
- Control Box: Easily accessible beneath the back rest panel for efficient servicing
- Quick Connect Bed Ends: Simple and quick to set-up upon bed installation. Easily

removed for access or maintenance

- Wall Buffer: Length adjustable to suit room configuration
- Cable Conduit: Allows other cables like air mattress power cords to be safely tucked away, preventing trip hazards and cord damage
- EasyMove Transport System: Facilitates safe no-lift transportation and movement of beds, as well as safe and easy access underneath the bed platform for routine inspection and maintenance
- Width Extension Kit: Ability to increase a single to a king single, ideal for resident changes in a facility
- High Performance Twin Castors x 8: Head end total locking (x2), foot end central locking (x4)

## Warranty Information

- 10 Years Frame
- 3 Years Electrics
- 1 Year Handset & Cables

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
NZ-BEB048600	Aidacare FLX Floorline Bed – Single	2130 - 2330 mm	964 mm	100 - 750 mm	275 kg	105 kg
NZ-BEB048610	Aidacare FLX Floorline Bed – King Single	2130 - 2330 mm	1124 mm	100 - 750 mm	275 kg	115 kg