

eCart Economy Carts



- Lightweight composite aluminium/plastic panels, reduces cart noise while in motion and panels can be replaced economically if damaged
- Fixed height push handle
- Recessed storage compartment
- Keyed mechanical lock
- Recessed ABS formed plastic top
- Slide out shelf fitted with runners for smooth, stable operation
- Raised drawer handles designed for grip and easy cleaning
- Drawers fitted with soft close runners and damper - prevents drawers from opening while in motion
- Epoxy powder coated frame
- Diameter 125 anti-static total lock castors
- Directional lock castor

Warranty Information

- 2 Years Electronics
- 5 Years Frame

Code	Display Name	Total Length	Total Width	Total Height
NZ-WAT901110	Economy Cart- 3 Drawer + Dividers (12) + Fixed Push Handle	520 mm	760 mm	1075 mm