



Wheelchair Accessories





Anti-tip Bars

Ensures chair does not tilt over when traveling across inclines and uneven surfaces. Comes in pairs.

Fits Socialite • Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward • Rehab RX2 • Rehab RX2 JNR



Attendant Brake Kits

Transit conversion Kit includes 12" wheels and attendant brakes, or attendant brakes kit without wheels.

Fits Evoke 2



Safety Lap Belt

Suitable for supporting users with low tone or weakness. Metal push button for secure fit.

Fits Vida • Vida X • Vida MG • Socialite • Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward • Rehab RX2 • Rehab RX2 JNR



Contoured Padded Overlay

Increased comfort, ideal for users in their wheelchair for longer periods of time. Breathable washable cover.

Fits Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward



IV Pole

Facilitates intravenous infusion.

Fits Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward • Rehab RX2 • Rehab RX2 JNR



Extension Lever

Extends brake lever for user accessibility.

Fits Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward



Calf Strap

Provides extra support for user's calves.

Fits Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward • Rehab RX2 • Rehab RX2 JNR



Base Frame Extension Kit

Extends stump support to cater for taller users.

Fits Rehab RX2 • Rehab RX2 JNR



Headrest Assembly

Deluxe angle adjustable head support. Helps to promote user posture and comfort.


Fits Evoke 2



Gutter Armrest Pad

Provides extra support for users with weak or unsupported arms.*

Fits Evoke 2



Adjustment Kit

Allows for easy, tool-free adjustment of a wheelchair.

Fits Rehab RX2 • Rehab RX2 JNR



Meal Tray

Provides sturdy surface for meal time.


Fits Rehab RX2



Mobility Accessory Bags

Keep personal items close and secure when you're on the move.

Fits Socialite • Vida • Vida X • Vida MG • Evoke 2 • Lite • Transit 3 • Assist 3



Leopard WAA696500



Floral WAA696510



Plaid WAA696520



Blue WAA696530



Drum Brake Kit

This upgrade kit can be retrofitted, providing the attendant with more control on sloped surfaces.


Fits Evoke 2 • Evoke 2 JNR • Rehab RX2 • Rehab RX2 JNR



Oxygen Bottle Holder

Holds oxygen bottle in place. Designed to hold C size bottles.

Fits Evoke 2 • Evoke 2 JNR • Rehab RX2 • Rehab RX2 JNR • Transit 3 • Assist 3



Stump Support

For amputee support.

(A) Fits Evoke 2 • Evoke 2 JNR • Rehab RX2 • Rehab RX2 JNR • Transit 2 • Assist 2


(B) Fits Transit 3 • Assist 3



Elevating Legrest

Padded legrests for improved comfort.

Fits Evoke 2 • Evoke 2 JNR • Transit 3 • Assist 3 • Assist 3 Ward • Rehab RX2 • Rehab RX2 JNR



Depth & Height Adjustable Armrests

Provides multiple support positions for users arms.

Fits Evoke 2



Extended Back Canes

Extends back canes up to 50mm to facilitate taller carers. Comes in pairs.

Fits Evoke 2 • Evoke 2 JNR



Stroller Handle

Angle adjustable handle allows greater ease of use for a carer.

Fits Rehab RX2 • Rehab RX2 JNR



Lateral Chest Support

Facilitates postural positioning.

Fits Evoke 2 • Evoke 2 JNR • Rehab RX2

Aspire Rehab RX2 Wheelchair Accessories

Included as Standard

Compatible

Code	Description	Rehab RX2	Rehab RX2 JNR
Stump Supports			
WCA848252	Left & Right - 410 to 460mm	<input checked="" type="radio"/>	<input checked="" type="radio"/>
WCA848254	Left & Right - 510 to 610mm	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Lateral Chest Support Assembly			
WCA848450	Left	<input checked="" type="radio"/>	<input type="radio"/>
WCA848460	Right	<input checked="" type="radio"/>	<input type="radio"/>
Elevating Legrests			
WCA848100	Left - 300/320/350/400/410mm	<input checked="" type="radio"/>	<input checked="" type="radio"/>
WCA848110	Right - 300/320/350/400/410mm	<input checked="" type="radio"/>	<input checked="" type="radio"/>
WCA848120	Left - 450/460/500/510mm	<input checked="" type="radio"/>	<input type="radio"/>
WCA848130	Right - 450/460/500/510mm	<input checked="" type="radio"/>	<input type="radio"/>
Base Frame Extension Kits			
WCA848530	50mm - Black	<input checked="" type="radio"/>	<input type="radio"/>
WCA848540	100mm - Black	<input checked="" type="radio"/>	<input type="radio"/>

Code	Description	Rehab RX2	Rehab RX2 JNR
Rear Wheels with Drum Brake Kits			
WCA848510	16"	<input checked="" type="radio"/>	<input checked="" type="radio"/>
WCA848520	20"	<input checked="" type="radio"/>	<input type="radio"/>
WCA848500	24"	<input type="radio"/>	<input type="radio"/>
Other			
WAA695600	Oxygen Bottle Bag	<input checked="" type="radio"/>	<input checked="" type="radio"/>
WCA848550	Tool-free Adjustment Kit	<input checked="" type="radio"/>	<input checked="" type="radio"/>
WCA848570	Meal Tray - 460mm to 510mm	<input checked="" type="radio"/>	<input type="radio"/>
WCA848270	Anti-tip Bars	<input type="radio"/>	<input type="radio"/>
WCA848280	IV Pole	<input checked="" type="radio"/>	<input checked="" type="radio"/>
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle	<input checked="" type="radio"/>	<input checked="" type="radio"/>
WCA848580	Calf Strap	<input checked="" type="radio"/>	<input checked="" type="radio"/>
WCA848590	Stroller Handle	<input checked="" type="radio"/>	<input checked="" type="radio"/>

*Note Side transfers not possible when gutter armrests are attached

Aspire Vida Range Wheelchair Accessories

Code	Description
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle
GSP047244	Anti Tippers

Aspire Socialite Wheelchair Accessories

Code	Description
WCA848270	Anti-tip Bars

Aspire Evoke 2, Transit 3 & Assist 3 Wheelchair Accessories

Code	Description	Evoke 2	Evoke 2 JNR	Transit 3	Assist 3
Contoured Padded Overlays					
WCA848325	300mm		●		
WCA848330	350mm	●			
WCA848340	400mm	●		●	●
WCA848350	450mm	●		●	●
WCA848360	500mm	●		●	●
WCA848370	570mm	●		●	●
WCA848380	610mm	●		●	●
Stump Supports					
WCA848220	Left - 300 to 410mm		●		
WCA848230	Right - 300 to 410mm		●		
WCA848252	Left & Right - 410 to 510mm	●		●	●
WCA848254	Left & Right - 510 to 610mm	●		●	●
WCA848260	Extension Bracket	●	●		
Elevating Legrests					
WCA848100	Left - 300 to 410mm	●	●		
WCA500080	Left - 400mm			●	●
WCA848110	Right - 300 to 410mm	●	●		
WCA500085	Right - 400mm			●	●
WCA848120	Left - 450 to 510mm	●			
WCA500070	Left - 450 to 610mm			●	●
WCA848130	Right - 450 to 510mm	●			
WCA500075	Right - 450 to 610mm			●	●
WCA848140	Left - 570/610mm	●			
WCA848150	Right - 570/610mm	●			
WCA848080	Left - Evoke 2 JNR		●		
WCA848090	Right - Evoke 2 JNR		●		
Lateral Chest Support Assembly					
WCA848430	Left	●	●		
WCA848440	Right	●	●		
Depth & Height Adjustable Armrests					
WCA848390	Left	●			
WCA848400	Right	●			

Code	Description	Evoke 2	Evoke 2 JNR	Transit 3	Assist 3
Gutter Armrest Pads*					
WCA848410	Left	●			
WCA848420	Right	●			
Rear Wheels with Drum Brake Kits					
WCA848510	16"		●		
WCA848520	20"		●		
WCA848500	24"	●			
Brake Kits					
WCA848592	Attendant Brakes Kit	●			
WCA848594	Transit Conversion Kit	●			
Brake Extension Levers					
WCA848310	Left	●	●	●	●
WCA848320	Right	●	●	●	●
Anti-tip Bars					
WCA848270	Anti-tip Bars	●	●		
WCA500050	Anti-tip Bars			●	●
WCA848272	Depth-adjustable Anti-tip Bars	●			
Other					
WCA848290	Oxygen Bottle Holder	●	●		
WCA848280	IV Pole	●	●		
WCA848285	IV Pole with Mounting Clamp			●	●
WCA848300	Safety Lap Belt	●	●		
WCA500055	Safety Lap Belt			●	●
WCA848470	Extended Back Cane 19"/21"	●	●		
WCA848480	Headrest Assembly 400 to 450mm	●	●		
WCA848481	Headrest Assembly 500 to 610mm	●			
WCA848580	Calf Strap	●	●		
WCA500065	Calf Strap			●	●
WCA500060	Tool Kit			●	●
AAL848590	Security Mod			●	●

*Side transfers not possible when gutter armrests are attached